 DE - Sicherheitshinweise - MEGASONEX.....	2
 EN - Safety Instructions - MEGASONEX.....	3
 IT - Istruzioni di Sicurezza - MEGASONEX.....	4
 FR - Consignes de Sécurité - MEGASONEX.....	6
 ES - Instrucciones de Seguridad - MEGASONEX	7
 SE - Säkerhetsanvisningar - MEGASONEX.....	9
 PL - Instrukcje Bezpieczeństwa - MEGASONEX.....	10
 NL - Veiligheidsinstructies - MEGASONEX.....	12

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EN - Safety Instructions - MEGASONEX

Declaration of Conformity

Goldspire Ltd. hereby declares that the product "NiMH Battery, No. 140141" complies with Regulation (EU) 2023/1542. The complete text of the EU Declaration of Conformity is available at: www.megasonex.com/Conformity.

As the user, you are responsible for the proper operation and handling of the Nickel-Metal Hydride (NiMH) batteries. Proper handling includes assembly, charging, usage, and selection of the application area. You are liable for any damages to the product itself or caused by it if the damages result from improper handling or usage errors. Therefore, please read the user manual carefully, as it contains important information and safety warnings.

⚠ WARNING! Nickel-Metal Hydride (NiMH) batteries can ignite/explode if used improperly, damaged, deeply discharged, or overcharged. Ensure proper use of the Nickel-Metal Hydride (NiMH) battery.

⚠ Electrical Hazards ⚠

As with any electrical device, to reduce the risk of electric shock:

- Keep the charging station and power adapter away from water or other liquids.
- Store the charging station or power adapter in a location where they cannot fall into or be pulled into water or another liquid. If this occurs, DO NOT touch the charging station but immediately unplug it. Allow the wet part to dry for several days before attempting to use it again. If water has entered the charging station, there may be corrosion of the internal components, which could cause the charging station to fail within months. Water damage to the charging station or power adapter is not covered by the warranty.

⚠ Warning Notes ⚠

To reduce the risk of burns, electric shocks, fire, or injury:

- Carefully supervise children or people with disabilities when using this product. Stop using it if you cannot handle it properly.

- ▶ Only use the device for its intended purpose, as described below. Do not use any accessories other than those provided by the manufacturer.
- ▶ Do not use the charging station, power adapter, handpiece, or brush heads if any component is loose or damaged. This product does not contain any repairable components.
- ▶ Keep all parts away from hot surfaces.
- ▶ Never use the power adapter if the power cord or plug is damaged: if the power cord is damaged, it must be replaced by the manufacturer or their customer service representative to avoid hazards.
- ▶ Unplug the power adapter before cleaning, and ensure it is completely dry before plugging it back in. Never fully immerse the charging station or power adapter in water or any other liquid.

Medical Warnings

- ▶ Consult your dentist before using this product if you have recently had oral surgery or gum treatment.
- ▶ Contact your dentist if heavy bleeding occurs during use or if mild bleeding persists after a week of use.
- ▶ If you have questions about using the product with a pacemaker or other implanted devices, consult your doctor or the device manufacturer before using the product. The MEGASONEX toothbrush has been tested and complies with safety standards for electromagnetic interference.
- ▶ Do not use the toothbrush if the brush head is broken or the bristles are bent, as damaged bristles could break off during brushing.
- ▶ This is an oral hygiene product designed for cleaning teeth and the tongue. Do not use it for any other purpose.
- ▶ This device should not be used by persons with reduced physical, sensory, or mental capabilities or by those lacking experience or knowledge unless supervised or instructed by a person responsible for their safety. The device has not been tested for the following groups and should not be used by them: pregnant women, diabetics, or patients with oral cancer.
- ▶ If you are unsure if there is any risk in using this product, always consult a qualified medical or dental professional.
- ▶ Children should only use the device under supervision to ensure it is not used as a toy.
- ▶ For hygiene reasons, wash the handpiece and any new brush heads with water and soap before the first use and do not share them with others.

Intended Purpose

The MEGASONEX toothbrush is designed for daily home use with any toothpaste of your choice. It can be used by individuals with braces, crowns, implants, veneers, and fillings. It is developed to combat gingivitis and plaque formation, which often lead to gum bleeding, gum recession, as well as tartar and cavities. Please note that the MEGASONEX toothbrush emits therapeutic ultrasound, but does not emit electromagnetic radiation, as is the case with cell phones or metal detectors.

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